

TECHNICAL DATA SHEET

Avalon FM30 Fetal Monitor with Intellispace Perinatal

Type : Avalon FM30 Fetal Monitor with Intellispace Perinatal
Parameter : Central Monitor for Fetal Monitor
COO :

Avalon FM30 Fetal Monitor with Intellispace Perinatal

US classifications : Pulsed doppler / 1 Mhz, 7 crystal
FHR coincidence detection capability : Yes

Output signals

Analog : No
Digital : Yes
Interated signal simulation / demo mode : Yes

Display Feature

DIGITAL FHR : <=3
Number of displays : One 16.5 cm (6.5") TFT color with touchscreen
Digital UA : Yes

Remote event marker : Yes

Special Feature

: Auto Intepretation NST
ISP Integration
Central Monitor for Fetal Monitor

Accessories

Toco MP Transducer 1
Two US Transducers 1
One Patient Module 1
Green Grid Paper 1
60mm Washable Belt 1
Toco Belt Button 1
TOCO Belt Clip 1