

## TECHNICAL DATA SHEET

### Philips Avalon Fetal Monitor FM 30 Cardiotocograph

Product	:	Avalon FM30 CTG
COO	:	Germany
US classifications	:	Pulsed doppler / 1 Mhz, 7 crystal
FHR coincidence detection capability	:	Yes
NST Trace Interpretation	:	Yes
Maternal parameter	:	SpO2 and NIBP
Output signals		
Analog	:	No
Digital	:	Yes
Interated signal simulation / demo mode	:	Yes
Display Feature		
DIGITAL FHR	:	<=3
Number of displays	:	One 16.5 cm (6.5") TFT color with touchscreen
Digital UA	:	Yes
Remote event marker	:	Yes
Accessories		
Toco MP Transducer	:	1
Two US Transducers	:	1
One Patient Module	:	1
Green Grid Paper	:	1
60mm Washable Belt	:	1
Toco Belt Button	:	1
TOCO Belt Clip	:	1
Battery	:	1
Trolley	:	1
SpO2 sensor	:	1
NIBP cuff	:	1