

TECHNICAL DATA SHEET Philips Defibrillator/AED Heartstart HS1

Type	:	Philips HeartStart Defibrillator HS1
COO	:	USA
Physical Approx	:	
	Size	Size : 7.2 cm H x 19 cm D x 21 cm W (2.8" H x 7.4" D x 8.3" W)
	Weight	With battery and pads cartridge: 1.5 kg (3.3 lbs.) Without battery or pads cartridge: 1 kg (2.4 lbs.)
Waveform	:	Truncated Exponential Biphasic; waveform parameters adjusted as a function of each patient's impedance.
Technology	:	Smart Biphasic and Smart Analysis
CPR guidance	:	Instructions for infants and children under 25 kg or 55 lbs. or 0-8 years old, and adults and children over 25 kg or 55 lbs or greater than 8 years old.
Battery	:	
	Type	9 Volt DC, 4.2 Ah, composed of disposable long-life lithium manganese dioxide primary cells.
	Capacity	Minimum 200 shocks or 4 hours of operating time.
	Standby life	Typically, 4 years when battery is installed and when stored and maintained according to directions provided in this document.
Accessories	:	
Battery		Include
Pads		Include
Case		Include