



TECHNICAL DATA SHEET

Philips Efficia Patient Monitor CM100

Type	: Efficia CM100
Parameter	: NBP, SpO2, HR, RR and Temperature
COO	: China
FDA approved	: Yes (K151812)
CE marks	: Yes
Screen size (apprx)	: apprx 10.1" (1280x800 pixels)
Display Rate	: TFT LED Display
Screen	: Touchscreen
Weight (apprx.)	: Approx. 3.3 kg
Waves	: Yes
Battery (apprx)	: Rechargeable Li-Ion battery (Maximum 2.5 hour continuous monitoring)
Trend review	: Yes
Full disclosure	: Yes
Standard measurement	: ECG/ST/QT/RESP/ SPO2
ECG	: EASI, QT Analysis and ST
RESP	: Impedancemetry lead II or lead I
SpO2	: FAST SpO2
NIBP	: Oscillometric with sequence mode and venous puncture
Temperature	: Yes
Accessories (for adult by default)	
SPO2 sensor	1
ECG cable	1
Temperature probe	1
NIBP cuff	1