



TECHNICAL DATA SHEET

Philips Efficia Patient Monitor CM150

Type	: Efficia CM150
Parameter	: NBP, SpO2, HR, RR, Temperature, EtCO2 and dual IBP
COO	: China
FDA approved	: Yes (K151812)
CE marks	: Yes
Screen size (apprx)	: apprx 15.6" (1366x768 pixels)
Display Rate	: TFT LED Display
Screen	: Touchscreen
Weight (apprx.)	: Approx. 5 kg
Waves	: Yes
Battery (apprx)	: Rechargeable Li-Ion battery (Maximum 2.5 hour continuous monitoring)
Trend review	: Yes
Full disclosure	: Yes
Standard measurement	: ECG/ST/QT/RESP/ SPO2
ECG	: EASI, QT Analysis and ST - 12 Lead ECG
RESP	: Impedancemetry lead II or lead I
SpO2	: FAST SpO2
NIBP	: Oscillometric with sequence mode and venous puncture
Temperature	: Yes
IBP	: Yes, dual IBP
EtCO2	: Yes
Accessories (for adult by default)	
SPO2 sensor	1
ECG cable	1
Temperature probe	1
NIBP cuff	1