



TECHNICAL DATA SHEET

Philips Goldway Patient Monitor G40E

Type	: Goldway G40E
Parameter	: NBP, SpO2, HR, RR and Temperature
COO	: China (CE IEC 4th/CFDA class III)
Screen size (apprx)	: 12.1"
Fan	: Fanless
Screen	: Touchscreen
Weight (apprx.)	: Approx. 4.2 kg
Waves	: Yes
Battery (apprx)	: approx 2.5 hours
Trend review	: Yes
Standard measurement	: ECG/ST/QT/RESP/ SPO2
ECG	: EASI, QT Analysis and ST
RESP	: Impedancemetry lead II or lead I
SpO2	: FAST SpO2
NIBP	: Oscillometric with sequence mode and venous puncture
Temperature	: Yes
Accessories (for adult by default)	
SPO2 sensor	1
ECG cable	1
Temperature probe	1
NIBP cuff	1