

TECHNICAL DATA SHEET

Philips Heartstart FRx Defibrillator (AED)

Type	:	Philips HeartStart FRx Defibrillator (AED)
COO	:	USA
Physical Approx	:	
	Size	18 x 6 x 22 cm (7.1 x 2.4 x 8.8 in) H x D x W
	Weight	Approximately 1.6 kg (3.5 lb) with battery and pads installed
		Biphasic truncated exponential; waveform
Waveform	:	parameters are automatically adjusted as a function of patient defibrillation impedance
		Verbal instructions for adult and infant/child
CPR guidance	:	CPR provides instructions and audio cues for the appropriate number, rate and depth of chest compressions, as well as for each breath
Battery	:	
	Type	9 Volt DC, 4.2 Ah, composed of disposable long-life lithium manganese dioxide primary cells.
	Capacity	Minimum 200 shocks or 4 hours of operating time.
	Standby life	Typically, 4 years when battery is installed and when stored and maintained according to directions provided in this document.
Accessories	:	
Battery		Include
Pads		Include
Case		Include