



TECHNICAL DATA SHEET

Philips Holter Digitrak XT

Type	: Philips Holter Digitrak XT
COO	: USA
Weight without battery apprx	: 62g
Width, height, depth	: apprx W 3.6in (91.44mm), H 2.2in (55.88mm), D 0.75in (19.05mm)
Display size	: 44mm x 22mm
Full disclosure	: Yes
Retrospective scanning	: Yes
Holter recorder	: 96 hours (1) and 24 hours (1)
EASI 12-lead ECG	: Yes
Pacer analysis	: Yes
Patient cable	: 1 Short, 1 Medium and 1 Long Cable
Battery life	: Single AAA alkaline battery for recording up to 7 days
Accessories	: Pouch, dual docking station, Belt Clip and electrode
Holter software application	: include with CD