

TECHNICAL DATA SHEET

Intellivue MX100 Patient Monitor

Type	: Intellivue MX100
Parameter	: NBP, SpO2, HR, RR, Temperature
COO	: Germany
Screen size (apprx)	: 6,1 inch
Size (W x H X D) in mm (apprx)	: 249 x 97 x 111 mm
Display Rate	: Active matrix color LCD display with capacitive multi-touch screen
Screen	: Antimicrobial Corning® Gorilla® Glass
Weight (apprx.)	: Approx. 1.4 kg
Waves	: Yes
Battery (apprx)	: 3-4 hrs
Smart key	: Yes and configurable
Trend review	: Yes
Standard measurement	: ECG/ST/QT/RESP/ SPO2
ECG	: EASI, QT Analysis and ST
RESP	: Impedancemetry lead II or lead I
SpO2	: FAST SpO2
NIBP	: Oscillometric with sequence mode and venous puncture
Temperature	: Yes
Accessories (for adult by default)	
SPO2 sensor	1
ECG cable	1
Temperature probe	1
NIBP cuff	1