

## TECHNICAL DATA SHEET

### Intellivue MX400 Advance

Type	:	Intellivue MX400
Parameter	:	NBP, SpO2, HR, RR, Temperature, IBP, EtCO2
COO	:	Germany
Screen size(approx.)	:	9" inch ; 800x480 WVGA
Size (W x H x D) in mm(apprx.)	:	277 x 288 x 186
Touch screen	:	Yes
Battery Life (approx.)	:	5 hours
Ambient light settings	:	Yes
Waves	:	4
Integrated recorder	:	Yes
Alarm indicator	:	Visual LED on the front
Standard measurement extension	:	Compatible with all Philips Multi Module Measurement Server
Multi Module Measurement Server		
Standard module	:	Yes
Standard measurement	:	ECG/ST/QT/RESP/ SPO2/IBP/EtCO2
ECG	:	EASI, QT Analysis and ST
RESP	:	Impedancemetry lead II or lead I
SpO2	:	FAST SpO2
NIBP	:	Oscillometric with sequence mode and venous puncture
Temperature	:	Yes
IBP	:	Yes
EtCO2	:	Yes
Accessories (for adult by default)		
SPO2 sensor	:	1
ECG cable	:	1
Temperature probe	:	1
NIBP cuff	:	1
EtCO2 Accessories	:	1