

## TECHNICAL DATA SHEET

### Intellivue MX550 Autocharting Patient Monitor

Type : Intellivue MX550  
 Parameter : NBP, SpO2, HR, RR, Temp, IBP, EtCO2, CO, BIS,  
 : Autocharting  
 COO : Germany

#### PHILIPS Intellivue MX550

Screen size (apprx) : 22 inch widescreen  
 Size (W x H X D) in mm (apprx) : (21.4 x 15.3 x 8.5 in)  
 Display Rate : 547 mm active matrix color LCD (TFT)  
 Screen : Touchscreen, 1920 x 1080 (Full HD)  
 Weight (apprx.) : 11 kg (24.3 lb)  
 Waves : > 8 waveform  
 Smart key : Yes  
 Trend review : 50 or 100 numerics @ 12 seconds, 1 minute, 5-minute  
 : resolution

Standard measurement : ECG/NIBP/RESP/ SPO2/TEMP/IBP/EtCO2/CO/BIS  
 ECG : ECG 5 lead, include EASI Method  
 RESP : Impedancemetry lead II or lead I  
 SpO2 : FAST SpO2  
 NIBP : Oscillometric with sequence mode and venous  
 : puncture  
 Temperature : Yes  
 IBP : Yes  
 EtCO2 : Yes  
 CO : Yes  
 BIS : Yes  
 Special Feature : AD 85  
 Special Feature : Autocharting  
 : Digital Flowsheet

Accessories  
 (for adult by default)

SPO2 sensor	1
ECG cable	1
Temperature probe	1
NIBP cuff	1
Mainstream sensor	1

#### AUTOCHARTING SYSTEM INTELLISPACE CRITICAL CARE & ANESTHESIA

Digital Flowsheet : Yes, Philips can manage based on hospital request

Calculation Engine : Yes, Philips can manage based on hospital request  
Clinical Advisories : Yes, Philips can manage based on hospital request  
Device Interfacing : Yes, for Patient Monitor  
Interoperability : Yes  
Notes & Forms : Yes  
Microbiology & pathology Results: Yes , if LIS in hospital available  
Medication and infusion manager: Yes, if CPOE in hospital available  
Self Service Clinical Intelligence : Optional