

## TECHNICAL DATA SHEET

### PIC iX Central Monitor (16 Beds)

Type : PIC iX Central Monitor  
Parameter : 16 bed license  
COO : USA

#### PHILIPS Information Center PIC iX

##### Hardware

Display : 19 inch (display size and type of monitor may vary)  
Single or dual display  
Hardware platform : HP/Compaq (PHILIPS standard)  
Switch : 1

##### Display information

Patient sector display : 12 - 16 (depends on request & capability)  
Graphical trends : Yes, split trends and ST-MAP within ST Sector  
Tabular trends : Yes  
Trends upload from monitor : Yes

Arrhythmia annalysis : 25 enhanced arrhythmia alarms  
ST analysis : Yes  
Full disclosure : Yes

##### Alarms

Type : High and low for all monitored parameters and arrhythmias  
Priority/ levels : Cyan \*, yellow \*\*, red \*\*\*, technical  
Silence active alarms : Yes  
All alarms suspend : Yes

##### Reports

Data type : Trend, alarm, event, and wave