

TECHNICAL DATA SHEET

PIC iX Central Monitor (8 Beds)

Type : PIC iX Central Monitor
Parameter : 8 bed license
COO : USA

PHILIPS Information Center PIC iX

Hardware

Display : 19 inch (display size and type of monitor may vary)
Single or dual display
Hardware platform : HP/Compaq (PHILIPS standard)
Switch : 1

Display information

Patient sector display : 12 - 16 (depends on request & capability)
Graphical trends : Yes, split trends and ST-MAP within ST Sector
Tabular trends : Yes
Trends upload from monitor : Yes

Arrhythmia annalysis : 25 enhanced arrhythmia alarms
ST analysis : Yes
Full disclosure : Yes

Alarms

Type : High and low for all monitored parameters and arrhythmias
Priority/ levels : Cyan *, yellow **, red ***, technical
Silence active alarms : Yes
All alarms suspend : Yes

Reports

Data type : Trend, alarm, event, and wave